

HI 3220 • HI 3221 • HI 3222

Laboratory pH Meters*with Backlit, Graphic LCD***Calibration Check™ pH/mV/ISE/Temperature Bench Meters**

HANNA's HI 3000 series of professional bench meters are designed to provide high accuracy and ease of use both in the laboratory as well as in harsh industrial conditions.

These meters feature HANNA's exclusive Calibration Check™, a diagnostics system that ensures accurate readings every time. By alerting users of potential problems during the calibration process, the Calibration Check™ system eliminates erroneous readings due to dirty or faulty pH electrodes or contaminated buffer solutions. Throughout the calibration process, users are guided step-by-step by the on-screen tutorial. After calibration, a probe condition indicator is displayed on the LCD informing the user of the overall electrode status.

These instruments can be calibrated at up to 5 points with 7 standard buffers and an additional 5 custom buffers. HI 3222, HI 3221 and HI 3220 also feature a user selectable "Outside Calibration" warning as well as a user selectable "Calibration Time-out" warning.

HI 3222, HI 3221 and HI 3220 feature a powerful interactive user support menu system that assists you before, during and after measurement. On-screen tutorials guide users through set-up, calibration and measurement while context sensitive help of any screen is available at a push of a button. The help screen includes

language specific assistance for menu parameters, calibration, log, contact information and accessories for your instrument.

HI 3000 series measures ORP with a resolution of 0.1 mV and also offers an extended temperature range from -20 to 120°C (-4 to 248°F) using HI 7662-T interchangeable probes.

HI 3221 and HI 3222 feature ISE measurement modes: HI 3221 displays ISE as ppm while the HI 3222 provides a choice of units (ppb, ppm, molarity, weight/volume %). The electrode type and unit selection capability of the HI 3222 and the ion charge selection capability of the HI 3221 make these instruments extremely versatile for a wide range of ISE measurements.

HI 3222, HI 3221 and HI 3220 feature three types of logging modes (for pH, ISE and mV): stability logging, interval logging and log-on-demand. Stability logging stores measurements once the electrode has stabilized in the solution. The user can set sensitivity between fast, medium and accurate. Interval logging stores data at user defined set time intervals and log-on-demand stores data at a touch of a button. Stored data includes the parameter, test results, offset and slope, date and time. The HI 3000 series can be connected to a PC via USB. HI 92000 Windows® compatible software (optional) aids in data management.

- Single (HI 3220 and HI 3221) or dual (HI 3222) input channels
- Backlit, Graphic LCD
- Exclusive Calibration Check™
- Prompts and messages on the graphic LCD
- Stability, interval and Log-on-demand logging
- 5 point calibration with 7 standard and 5 custom buffers
- Relative mV measurements

- Up to 400 sample memory capacity
- Auto HOLD to freeze readings on the display
- User selectable language
- GLP to view last calibration data for pH, Rel mV or ISE
- User selectable "Outside Calibration Range" warning
- User-selectable "Calibration Timeout" to remind when a new calibration is necessary
- PC interface via USB

HI 3220 • HI 3221 • HI 3222

Laboratory pH Meters

with Backlit, Graphic LCD

5

ISE/pH




Measurement

The HI 322x series large backlit graphic display shows multiple messages along with the current measurement readings simultaneously.

Screen specific commands can be easily administered via soft keys located on the front panel.


Help menus, set-up screens and more can be accessed through the button panel on the right.

Measurement Features



pH Measurement

pH readings are displayed along with temperature, calibration points, electrode condition, selected temperature compensation.




ISE Measurement

An ISE electrode can be attached to the HI 3221 and HI 3222.




ORP Measurement

Just connect an optional ORP electrode to perform ORP measurements. Measurements are displayed with 0.1 mV resolution.



Relative mV Measurement

These instruments can also display Relative mV. When reading, the Relative mV will be displayed along with the Absolute mV value and the current temperature.




Temperature Measurement

Use the HI 7662-T temperature probe to automatically compensate for pH measurements. Measurements can be displayed in °C or °F.

Laboratory pH Meters


with Backlit, Graphic LCD

HANNA's Exclusive Calibration Check™ Features




Calibration

Detailed messages throughout the easy pH, ORP, Relative mV and ISE calibration process.





Electrode Condition

Electrode status is indicated on the display after calibration.




Wrong Buffer Warning

The "Wrong Buffer" warning is displayed when calibration cannot be confirmed.




Wrong Slope Warning

The "Wrong Slope" warning message is displayed if the slope is out of the accepted range.




Clean Electrode Warning

The "Clean Electrode" warning message is displayed to advise the user that the pH electrode may need to be cleaned.



Electrode Dirty/Broken and Buffer Contaminated Warnings


These warnings are displayed when the offset of the electrode is not in the accepted range and calibration cannot be confirmed. This could be the result of a broken electrode, or an electrode which needs to be cleaned. Also, the buffer solution quality could have been compromised and may need to be changed.



GLP


Last calibration data is stored automatically after a successful calibration. Calibration data can be easily accessed and includes calibration buffer, offset, slope and electrode condition information. Previous calibration points are shown in reverse text.

Logging




Log and View Results

This feature allows users to log pH, Rel mV or ISE measurements for later retrieval. All logged data can be transferred to a PC through USB.



Set-up



Set-up

General parameters can be easily modified through the set-up menus.

SPECIFICATIONS

HI 3220

HI 3221

HI 3222

Range	pH	-2.0 to 20.0; -2.00 to 20.00; -2.000 to 20.000 pH		
	mV	±2000 mV		
	ISE	—	From 1.00 x 10 ⁻³ to 1.00 x 10 ⁵ ppm	From 1.00 x 10 ⁻⁷ to 9.99 x 10 ¹⁰ concentration (choice of units)
	Temperature	-20.0 to 120.0 °C (-4.0 to 248.0°F)		
Resolution	pH	0.1; 0.01; 0.001 pH		
	mV	0.1 mV		
	ISE	—	3 digits 0.01; 0.1; 1; 10 concentration	
	Temperature	0.1°C (0.1°F)		
Accuracy	pH	±0.01; ±0.002 pH		
	mV	±0.2 mV		
	ISE	—	±0.5% of reading (monovalent ions), ±1% of reading (divalent ions)	
	Temperature	±0.2°C (±0.4°F) (excluding probe error)		
Calibration	pH	Up to 5 point calibration, 7 standard buffers available (1.68, 4.01, 6.86, 7.01, 9.18, 10.01, 12.45) + 5 custom buffers		
	ISE	—	Up to 2 point calibration, 6 standard solutions (0.1, 1, 10, 100, 1000, 10000 ppm)	Up to 5 point calibration, 6 standard solutions (in units selected)
	Slope	From 80 to 110%		
Temperature Compensation (pH)		Manual or automatic from -20.0 to 120.0°C (-4.0 to 248.0°F)		
Input Channels		1	1	2
Probes		HI 1131B pH with BNC connector and HI 7662-T stainless steel temperature probe (included)		
Logging		Log-on-demand 200 samples	Log-on-demand 300 samples	Log-on-demand 400 samples
Lot Logging		5, 10, 30 sec 1, 2, 5, 10, 15, 30, 60, 120, 180 min (max 600 samples)		
PC Connectivity		Opto-isolated USB (with HI 92000 software)		
Input Impedance		10 ¹² Ohms		
Power Supply		12 VDC		
Environment		0-50°C (32 to 122°F) Max. RH 55% non-condensing		
Dimensions		235 x 207 x 110 mm (9.2 x 8.14 x 4.33")		
Weight		1.8 kg (4.1 lbs.)		

ORDERING INFORMATION

HI 3220-01 (115V) and **HI 3220-02** (230V), **HI 3221-01** (115V) and **HI 3221-02** (230V), **HI 3222-01** (115V) and **HI 3222-02** (230V) are supplied with HI 1131B pH electrode with BNC connector and 1 m cable, HI 7662-T stainless steel temperature probe, pH 4 and 7 buffer solutions, cleaning solution, electrolyte solution, electrode holder, 12 Vdc power adapter and instructions.

ELECTRODES

All electrodes part numbers ending in "B" are supplied with a BNC connector & 1 m (3.3') cable, as shown below:

- HI 1043B** Use: *strong acid/alkalis*; Glass-body, double junction, refillable, combination pH electrode
- HI 1053B** Use: *emulsions*; Glass-body, triple ceramic, refillable, combination pH electrode
- HI 1083B** Use: *biotechnology*; Glass-body, open junction, refillable, combination pH electrode

HI 1131B

Use: *general purpose*; Glass-body, single junction, refillable, combination pH electrode

HI 3230B

Use: *general purpose*; Plastic-body, gel-filled, combination ORP electrode

platinum

HI 7662-T

Stainless steel temperature probe

SOLUTIONS

HI 5004

pH 4.01 buffer solution, 500 mL

HI 5007

pH 7.01 buffer solution, 500 mL

HI 5010

pH 10.01 buffer solution, 500 mL

HI 54710

pH 4.01, pH 7.01 and pH 10.01 buffer solution, 500 mL ea.

HI 7020L

ORP test solution, 200-275 mV, 500 mL

HI 7021L

ORP test solution, 240 mV, 500 mL

HI 7022L

ORP test solution, 470 mV, 500 mL

HI 7091L

Reducing Pretreatment ORP solution, 500 mL

HI 7092L

Oxidizing Pretreatment ORP solution, 500 mL

HI 7071

3.5M KCl + AgCl Electrolyte, (4) 30 mL, for single junction electrodes

HI 7082

3.5M KCl Electrolyte, (4) 30 mL, for double junction electrodes

HI 7061L

Electrode cleaning solution, 500 mL

HI 70300L

Electrode storage solution, 500 mL

HI 50XX buffers are technical buffers with ±.01 accuracy and are provided with certificate

ACCESSORIES

HI 740157 Plastic refilling pipette (20)

HI 76405 Electrode holder

HI 92000 Windows® compatible software

HI 920013 USB cable for PC connection